

NittoPhase[®]HL Solid Support

High Loaded Polymeric Solid Supports for Oligonucleotide Synthesis



Mohammad Ahmadian, Tatsuya Konishi, Kenjiro Mori, Adrian Dukanovic, John Stevenson, Paula Gunawan, Yuki Kinjo, Eri Maeta, Liping Wang

Introduction

Increasing interest in the development of oligonucleotide therapeutic molecules, especially RNAi, has led to a growing need to optimize synthesis. Materials with higher performance to price ratio will allow drug development companies to meet their API cost targets. Solid support is one component that has come under close scrutiny in this regard. Controlled Pore Glass (CPG) supports have a limited loading capacity of around 80-90 $\mu\text{mol/g}$, while polymeric supports, capable of achieving significantly higher loading, have historically failed to show satisfactory and reliable performance at loadings greater than 200 $\mu\text{mol/g}$ for DNA and 150 $\mu\text{mol/g}$ for RNA synthesis. To address this issue, we have developed NittoPhase[®]HL, a high performance solid support with optimized loading capacities of up to 250 $\mu\text{mol/g}$ and 400 $\mu\text{mol/g}$ for RNA & DNA oligonucleotide syntheses, respectively. Since its development in 2010, NittoPhase[®]HL has shown superb performance in both synthesis yields and purities as compared to other solid supports in side-by-side comparisons. Regardless whether it is a 1-2 mmol lab scale or clinical trial material production batch (at 220, 550, 700 and 900 mmol scales), NittoPhase[®]HL has always produced highly favorable results in DNA synthesis. In these studies, NittoPhase[®]HL loaded with UnyLinker[™] at 350 $\mu\text{mol/g}$ was utilized for synthesis of 20-mer MOE Gappers and achieved per micromole yields ($\text{OD}_{260}/\mu\text{mol}$) higher or equivalent to alternative supports. Meanwhile, NittoPhase[®]HL at loading of 250 $\mu\text{mol/g}$ has also been successfully utilized in syntheses of highly modified siRNAs in various scales (ranging from 65 μmol to 65 mmol) to provide prominent crude yields and purities. These increased yields and purities, together with the elevated loading capacities, undoubtedly empower NittoPhase[®]HL to offer superior performance at lower unit cost for oligonucleotide synthesis.

DNA Synthesis Results*

Sequence	Synthesis Scale (mmol)	Support Type	Loading ($\mu\text{mol/g}$)	Crude Purity (%)	Crude Yield ($\text{OD}_{260}/\mu\text{mol}$)
5-10-5 MOE gapper	1	Support A	360	78.5	93
		NittoPhase [®] HL	431	78.6	93
5-10-5 MOE gapper	1	NittoPhase [®] HL	431	68.8	80
		NittoPhase [®] HL	350	80.4	101
5-10-5 MOE gapper	2	NittoPhase [®] HL	350	75.2	94
		NittoPhase [®] HL	400	71.4	90
5-10-5 MOE gapper	2	NittoPhase [®] HL	400	72.8	93
		Support A	360	82.8	100
3-10-3 MOE gapper	2	Support A	360	83.1	100
		NittoPhase [®] HL	324	82.8	101
3-10-3 MOE gapper	2	NittoPhase [®] HL	324	83.0	105
		Support A	316	81.0	108
5-10-5 MOE gapper	220	Support A	316	79.6	100
5-10-5 MOE gapper	220	Support A	325	79.6	100
5-10-5 MOE gapper	220	NittoPhase [®] HL	319	83.2	108
5-10-5 MOE gapper	550	NittoPhase [®] HL	319	86.3	111
5-10-5 MOE gapper	700	NittoPhase [®] HL	315	90.7	112
5-10-5 MOE gapper	900	NittoPhase [®] HL	315	83.2	105

* Data courtesy of Isis Pharmaceuticals, Inc.




Figure 1. IP-HPLC-UV Trace for NittoPhase[®]HL UnyLinker[™] DNA synthesis.

Standard RNA Synthesis Results

Sequence	Synthesis Scale	Support Type	Loading ($\mu\text{mol/g}$)	Weight Yield ($\text{mg}/\mu\text{mol}$)	Crude Purity (%)	Crude yield ($\text{OD}_{260}/\mu\text{mol}$)
5-5-5-1-mer	200 μmol	Support A	308	5.06	71.2	141
5-5-5-1-mer	175 μmol	NittoPhase [®] HL	251	6.22	75.4	175
3-5-5-6-2-mer	303 μmol	NittoPhase [®] HL	244	8.76	74.0	124
6-5-5-3-2-mer	567 μmol	NittoPhase [®] HL	207	9.49	84.8	112
4-2-7-6-2-mer	1 mmol	NittoPhase [®] HL	253	7.27	62.4	87
4-2-7-6-2-mer	4 mmol	NittoPhase [®] HL	253	8.56	77.4	120




Figure 2. IP-HPLC-UV Trace for NittoPhase[®]HL T250 RNA synthesis.




Figure 3. Custom loaded solid support performance comparison.

Highly Modified RNA Synthesis Results

Support Type	Loading ($\mu\text{mol/g}$)	Synthesis Scale	Strand	Crude Purity (%)	Crude Yield ($\text{OD}_{260}/\mu\text{mol}$)
Support C	30	65 μmol	Sense	61	122
			Antisense	58	78
NittoPhase [®] HL	250	325 μmol	Sense	74	128
			Antisense	68	81
NittoPhase [®] HL	250	2 mmol	Sense	79	145
			Antisense	84	142
NittoPhase [®] HL	250	16.4 mmol**	Sense	73	149
			Antisense	78	157
NittoPhase [®] HL	250	65 mmol**	Sense	76	151
			Antisense	79	136

** In collaboration with NITTO DENKO AVECIA Inc. & Quark Pharmaceuticals, Inc.

Polymer Support Properties

Product Properties	NittoPhase [®]	NittoPhase [®] HL
Functionality	Hydroxyl group	Hydroxyl group
# of Hydroxyl Group	420 $\mu\text{mol/g}$	550 $\mu\text{mol/g}$
Average Particle Size	90 μm	85 μm
Average Pore Size	55 nm	45 nm
Dry Volume	3.2 ml/g	2.7 ml/g
Optimal loading	200 $\mu\text{mol/g}$	200 – 400 $\mu\text{mol/g}$
Swelling Volume (in ACN)	3.85 ml/g	4.0 ml/g
Swelling Volume (in toluene)	6.4 ml/g	6.1 ml/g
Column Packing Recommendation (6.3 ml Column)	DNA	1.10 g (20-mer, 200 $\mu\text{mol/g}$)
	RNA	1.05 g (21-mer, 150 $\mu\text{mol/g}$)
		0.80 g (20-mer, 350 $\mu\text{mol/g}$)
		0.69 g (21-mer, 250 $\mu\text{mol/g}$)